



S H E L L E Y
M E D I C A L
I M A G I N G
T E C H N O L O G I E S

Sales Office

London, Ontario, Canada


Phone: 1 (519) 690-0874

Fax: 1 (519) 690-0875


Email: sales@simutec.com

Web: www.simutec.com

Silicone Introducer Guide



Inside the white silicone body is a guide that will guide your needle.



Simply insert the needle until you feel the plastic introducer and then insert the needle in the middle of it.

Once the needle is inserted you can put your catheter through. This will provide you with a leak free connection to the models.

ELSTRAT
Product Line