



S H E L L E Y
M E D I C A L
I M A G I N G
T E C H N O L O G I E S

S a l e s O f f i c e

London, Ontario, Canada

Phone: 1 (519) 690-0874

Fax: 1 (519) 690-0875

Email: sales@simutec.com

Web: www.simutec.com

Luer Introducer



The Luer Introducer lets you adapt all of your Luer devices to our Elastrat line of models. You can easily introduce your catheters, syringes, or Luer adapters without any leakage. Connects to 5/8 and 6/10 mm tubes.



ELASTRAT in vitro models are compatible with modern imaging modalities such as digital subtraction angiography, computed tomography and magnetic resonance imaging. Providing the use of an adequate circulating fluid, Doppler techniques can also be performed. The in vitro models transparency to light makes them suitable for video and photographic monitoring.

ELASTRAT
Product Line