




SHELLEY
MEDICAL
IMAGING
TECHNOLOGIES

Sales Office
London, Ontario, Canada
Phone: 1 (519) 690-0874
Fax: 1 (519) 690-0875
Email: sales@simutec.com
Web: www.simutec.com

Work-tray and Sawhorses



Our work-tray in polyethylene has been designed to be compatible with all vascular models displayed in this catalogue.



Adjustable sawhorses.
In height: from 77cm to 95cm (30" to 37.4").

In width: from 75cm to 120cm (29.5" to 47.2").